

Logo Placements 3014 Canvas+™ Shorts

Recommended logo placements Heat seal & Embriodery

- Right leg, 1 cm from stud on ruler pocket maximum size 120x120 mm
- 100 Right leg, centered 1 cm abovereflective seam, maximum size 120x120 mm
- 101 Right leg, 1 cm from stud on ruler pocket maximum size 80x80mm
- 18 Centred at left hip-pocket, maximum size 80x60 mm

Heat seal

18 Centred at right hip-pocket, maximum size 120x60 mm



Add
B after the placement and it will be below the (1B)
Y and it will be turned from bottom to top

Z after the placement and it will be turned from top to bottom

